

LifeViz[®] *Infinity* ∞



A UNIQUE,
PORTABLE 3D,
360° IMAGING SYSTEM
TO CONVINCe PATIENTS
TO ENGAGE IN
AESTHETIC PROCEDURES

1 SYSTEM

COUNTLESS POSSIBILITIES




A PORTABLE IMAGING
SYSTEM FOR THE FACE,
BODY AND BREAST.

Portability

The LifeViz® Infinity is a unique system composed of a 3D camera mounted on a tripod with a turntable allowing quick capture of a 360° 3D body image. In hand-held mode, the camera is capable of generating 3D images of face and breast to simulate procedures and photo-document before/after interventions, without blocking an entire exam room with cumbersome equipment.

CAPTURE HIGH QUALITY 3D IMAGES



Accurate & Reproducible

Simply click, align the dual beams to the proper distance and take pictures of face, body or breast areas without a repositioning device. The dual flash ensures consistent lighting over multiple visits and under varying conditions.

Quick & Easy

In under three minutes, view and assess your patient needs for face, body or breast. The LifeViz® system has been developed closely with practicing physicians to ensure a smooth fit into office workflow and with staff.



DIFFERENTIATE
YOUR PRACTICE



IMPROVE PATIENT
COMMUNICATION



360°
SHOW
RESULTS



CONVINCE
& REASSURE




INCREASE
PATIENT CONVERSIONS

1 ALL
IN
ONE



www.quantificare.com



ANGLES

MEASUREMENTS

719 mm

759 mm

811 mm

By courtesy of Dr Paolo Bonan (Italy)

PRE-OP

Treatment planning

QUICKLY OBTAIN MEASUREMENTS

Get baseline measures with one-click and recommend suitable procedures.

EASE TREATMENT SELECTION

Use real-time simulation to align with patient expectations tailored to their anatomy.

Make more effective use of contouring devices to anticipate volume changes.

For breast augmentation, select the ideally suited breast implants from your personalized catalogue and simulate realistic results.

ASSESS SKIN CONDITIONS TO PERSONALIZE TREATMENTS

Evaluate the patient's skin health using multivariate analysis for the following measurements:

- Wrinkles
- Pores
- Oiliness
- Evenness
- Red spots / Brown spots
- Red map / Brown map.

Generate 3D animations and a report of skin health parameters to guide the treatment approach and communicate suitable aesthetic options.



ACCURATELY COMPARE BEFORE & AFTER TREATMENT IMAGES


Choose among a variety of comparison options (slider, Silhouette, side by side) and export them in high quality images or animated videos.

Easily track results and compare circumference, volume and shape.


Improve patient follow-up with the 3D photo documentation software which includes pigmentation and vascularization evaluations.

IMPROVE COMMUNICATION AND PATIENT SATISFACTION


Use an objective assessment tool to check your artistry, reassure your patient about treatment outcomes and provide comprehensive recommendations for future interventions.



PRE-OP



POST-OP



VOLUME

0.57 mL

SLIDER

By courtesy of Dr Pierre Moullot (France)

POST-OP

Before/After comparison